2024-05-16 Signature:

Institute: AAC

## **Operating Instructions**

According to the Betriebssicherheitsverordnung



Faculty of Chemistry

Work place: AAC Labs

## EQUIPMENT

Interim gas storage (garage at refectory/S05 T01 B52) Argon bundle

#### DANGERS FOR HUMANS AND ENVIRNOMENT











- Gas cylinders generally contain gases in a compressed state (air, oxygen, nitrogen, hydrogen,...), or in a liquefied state (ammonia, chlorine, carbonic acid, propane, butane), or as a gas in/with a carrier medium (acethylene with acetone).
- With a liquid carrier medium, the gas cylinder must never be opened or operated lying down. It should be noted that large quantities of gas can be released even from small cylinders.
- Some gases are odourless and tasteless.
- Gases with a molecular weight > 32 collect near the ground.
- **Danger of suffocation**: If released spontaneously, gases displace air. This can lead to dangerously low oxygen concentrations in the breathing air without being noticed.
- In the case of oxygen, it can lead to spontaneous fires in combination with organic substances. High oxygen concentrations can also be dangerous for humans.
- Compressed gases can also have dangerous substance properties:
  - acutely toxic
  - special health hazards (**CMR** [carcinogenic, germ cell mutagenic, toxic for reproduction], **STOT** [specific target organ toxicity], sensitising...)
  - corrosive (or irritant)
  - flammable or pyrophoric
  - oxidising
- **Explosion hazard**: If flammable gases escape, explosive gas/air mixtures are often formed.
- Risk of poisoning if toxic gases escape and are inhaled.
- **Risk of burns**: If corrosive gases escape, severe damage is possible by inhalation or contact with eyes and skin.

#### PROTECTIVE MEASURES AND RULES OF CONDUCT





- Organisational: Carry out storage and retrieval only in accordance with the rules of use and the storage instruction.
- Personal:
- Wear protective goggles.
- Use protective leather gloves.
- If necessary, carry a respiratory protection mask (see instructions for use and storage).
- Technical:
- Use only suitable means of transport (cylinder trolley).









- Secure all gas cylinders with chains to prevent them from falling over.
- · Protect containers from heat.
- Leave the intermediate gas cylinder storage immediately if gases escape.
- Take injured persons out of the danger zone in selfprotection (observe substance properties and oxygen deficiency).
- If the gas begins to escape, try to close the main valve, paying attention to your own safety.
- Caution: Respiratory protection mask with suitable filter can only be used for pollutant concentrations up to 1 %!
- Ask endangered persons in the surrounding work areas to leave the danger area.
- In case of *contact with clothing*, remove affected clothing immediately and allow to air.
- In case of skin contact, treat according to the substance properties, usually wash off with plenty of water. If necessary, consult a doctor.
- In case of eye contact, treat according to the substance properties, as a rule rinse the eyes with water for at least 15 minutes. Immediately consult an ophthalmologist.
- In case of inhalation:
- High concentrations of chemically inert gases may cause suffocation. Symptoms: Loss of mobility and consciousness. The affected person does not notice the suffocation. Remove injured person to fresh air using selfcontained breathing apparatus. Keep warm and calm. In case of respiratory arrest, give artificial respiration. Call a doctor!
- Slightly water-soluble (corrosive or) irritant gases may cause pulmonary oedema. Use auxilosone spray immediately. Call a doctor!
- After inhalation of toxic or harmful gases, administer oxygen if possible. Call an emergency doctor!

# RESPONSE TO MALFUNCTIONS



Shut down device immediately, faults must be reported to:

Constanze Dietrich Tel.: 1834559
Dr. Sven Meckelmann Tel.: 1833295
Prof. Dr. Oliver Schmitz Tel.: 1833950

### BEHAVIOR IN CASE OF ACCIDENT/FIRST AID



- Keep calm.
- Call first responders.
- Emergency call: 0112
- Report accident.

- Maintenance only by authorized, competent persons.
- Switch off the device and disconnect it from the mains.
- Regular inspection of wearing parts.